

## soft air® cumulus wds

Stage III



Dynamic ulcer pressure relief system with large cells for prevention of ulcers and therapy support for medium to high risk/stage III








ADL® GmbH  
 An der Germania Brauerei 1  
 48159 Münster  
 Germany  
 Telefon: +49 (0) 251 149896-0  
 Telefax: +49 (0) 251 149896-99

ADL® GmbH  
 Produktionsstätte  
 und Servicecenter  
 Im Flürchen 1  
 99334 Amt Wachsenburg  
 Germany  
 Telefon: +49 (0) 36202 7519-0  
 Telefax: +49 (0) 36202 7519-12

E-Mail: [info@adl-gmbh.de](mailto:info@adl-gmbh.de)  
 Internet: [www.adl-gmbh.de](http://www.adl-gmbh.de)



### Technical details

	soft air® cumulus wds
Art. No.	110 400-wds-cumulus
Material	Mattress: PU · Pump: ABS housing · Cover: PU with 3 sided zipper
Warranty	2 years (wear and tear excluded)
Electrical values	230 V · 50 Hz/7 W
Fuse	500 mA/slow
Dimensions (inflated)	200 x 84,5 x 18,5 cm (79" x 33" x 7½")
Weight	Mattress: ca. 6,5 kg (14.5 lbs) · Pump: ca. 1,5 kg (3.5 lbs) · Cover: ca. 0,9 kg (2 lbs)
Cleaning	Thermal disinfection
Cover	Washable at up to 95 °C, dryable at up to 50 °C     
Recommended patient weight	30 – 140 kg (5 – 22 st)
Application area	Ulcer prevention and therapy support for medium to high risk/stage III